

Accu-Fit Back™ Reference/Measurement Guides

"BETWEEN THE RAILS" BACK		
Wheelchair Width	Accu-Fit Back Width	Chest Width (Min/Max)
12"	8.5"	6" - 9"
13"	9.5"	7" - 10"
14"	10.5"	8" - 11"
15"	11.5"	9" - 12"
16"	12.5"	10" - 13"
17"	13.5"	11" - 14"
18"	14.5"	12" - 15"
19"	15.5"	13" - 16"
20"	16.5"	14" - 17"
21"	17.5"	15" - 18"
22"	18.5"	16" - 19"

"FULL WIDTH" BACK		
Wheelchair Width	Accu-Fit Back Width	Chest Width (Min/Max)
12"	12"	9.5"-12.5"
13"	13"	10.5"-13.5"
14"	14"	11.5"-14.5"
15"	15"	12.5"-15.5"
16"	16"	13.5"-16.5"
17"	17"	14.5"-17.5"
18"	18"	15.5"-18.5"
19"	19"	16.5"-19.5"
20"	20"	17.5"-20.5"
21"	21"	18.5"-21.5"
22"	22"	19.5"-22.5"

Above measurements are based on standard foam thickness of 1.5 Inches. Chest width measurements are measured from center-to-center of Accu-Fit Laterals. All measurements/dimensions are approximate.

Updated 4/1/2013